The GeoAxon Kuduwave – A Validated Type 1 Clinical Audiometer for Comprehensive Diagnostic Assessment

Authors: Dr HL (Dirk) Koekemoer

Affiliation: GeoAxon

Keywords: Kuduwave, Type 1 Clinical Audiometer, Diagnostic Audiology, Immittance Audiometry, GeoAxon, Boothless Audiometry, Tympanometry, Acoustic Reflexes, IEC 60645-1, IEC 60645-5. **Date of Original Publication:** June 13, 2025

Date of Revision/Expansion: N/A

Version: 1.0 (Original publication)

Competing Interests Statement: Dr. HL (Dirk) Koekemoer has a competing interest to declare. He is a consultant for GeoAxon, the manufacturer of the Kuduwave audiometer. This work was also supported by funding from GeoAxon, as noted in the Funding Statement. The author has strived to provide an objective, evidence-based analysis.

Funding Statement: This work was supported by GeoAxon.

Abstract

This white paper aims to establish the GeoAxon Kuduwave (Prime, Pro, and Pro TMP models) as a Type 1 Clinical Audiometer, compliant with stringent international standards including IEC 60645-1 for pure-tone audiometry and IEC 60645-5 for immittance audiometry (Pro TMP model). It details the Kuduwave's extensive diagnostic capabilities, encompassing air and bone conduction (including extended high frequencies), speech audiometry, and, for the Pro TMP model, comprehensive immittance testing including tympanometry (226 Hz probe tone, bilateral simultaneous capability) and acoustic reflexes (ipsilateral, contralateral, and bilateral). The paper synthesizes evidence from GeoAxon's internal certification and testing documentation, alongside a robust body of peer-reviewed publications and clinical validation studies. It addresses the evolution from earlier eMoyo-manufactured Kuduwave models, which included Type 3 diagnostic audiometers, to the current advanced GeoAxon platform as a Type 1 clinical audiometer. The Kuduwave's performance, validated in both traditional sound booth environments and innovative boothless settings, underscores its accuracy, reliability, and versatility. This document provides compelling evidence for audiologists and ENT surgeons that the GeoAxon Kuduwave is a sophisticated, reliable, and comprehensive diagnostic tool for contemporary audiological practice.

1. Introduction: Advancing Diagnostic Audiology with the GeoAxon Kuduwave

The Modern Imperative for Comprehensive and Accessible Audiological Diagnostics:

The global burden of hearing loss is substantial and projected to increase, with estimates suggesting nearly 900 million people could be affected by 2050.¹ This escalating prevalence underscores the critical need for accurate, timely, and broadly accessible audiological diagnostic services.² Traditional audiological practices, however, face inherent challenges. The reliance on expensive, immobile sound-treated booths and the scarcity of specialist hearing healthcare professionals, particularly in low- to middle-income countries and remote regions, significantly limit service delivery and accessibility.⁵ These limitations necessitate innovative solutions that can overcome geographical and financial barriers without compromising diagnostic quality.

Overview of the GeoAxon Kuduwave as an Innovative Solution:

The Kuduwave audiometer, manufactured by GeoAxon Global (PTY) Ltd, emerges as a significant technological advancement designed to address these contemporary challenges in audiology. It is a portable, computer-based audiometric system engineered with core philosophies of boothless operation, comprehensive diagnostic capabilities, and tele-audiology enablement. The Kuduwave's design moves beyond the physical constraints of traditional testing environments. Its validation for use outside of sound-treated rooms fundamentally alters the paradigm of how and where high-quality audiological diagnostics can be performed. This capability has profound implications for enhancing accessibility in diverse settings, including remote or rural areas, mobile clinics, and even within space-constrained hospital environments, thereby potentially democratizing access to advanced hearing diagnostics.

Furthermore, the evolution of the Kuduwave platform from earlier models manufactured by eMoyo to the current GeoAxon systems is noteworthy. This progression indicates a focused commitment to advancing the technology to meet and exceed higher clinical standards. While older iterations had their own classifications, the current GeoAxon Kuduwave models have achieved more stringent certifications, reflecting a deliberate investment in engineering and compliance to deliver premier clinical-grade performance. Understanding this technological trajectory is crucial, as it highlights a dedication to quality and continuous improvement, addressing any potentially outdated perceptions based on earlier versions of the device.

Statement of Purpose:

This white paper aims to provide a comprehensive, evidence-based demonstration of the GeoAxon Kuduwave's classification as a Type 1 Clinical Audiometer. It will detail the full suite of diagnostic functionalities offered by the Kuduwave Prime, Pro, and Pro TMP models. The objective is to furnish audiologists and Ear, Nose, and Throat (ENT) surgeons with compelling data and analyses, assuring them of the Kuduwave's suitability and reliability for comprehensive audiological assessments in contemporary clinical practice.

2. Establishing the GeoAxon Kuduwave as a Type 1 Clinical Audiometer

The classification of audiometric equipment is governed by rigorous international standards that define the accuracy, range, and capabilities required for different levels of audiological assessment. A Type 1 classification represents the highest echelon, signifying suitability for comprehensive clinical and diagnostic evaluations. This section details the GeoAxon Kuduwave's adherence to these premier standards.

Defining International Standards for Clinical Audiometers and Immittance Instruments:

The International Electrotechnical Commission (IEC) sets forth key standards for audiometric equipment. IEC 60645-1, "Electroacoustics – Audiometric equipment – Part 1: Equipment for pure-tone audiometry," categorizes audiometers into Types 1 through 4. Type 1 audiometers are designated for clinical and diagnostic use, demanding the most stringent performance in terms of accuracy, frequency range, intensity range, and available test signals. Similarly, IEC 60645-5, "Electroacoustics – Audiometric equipment – Part 5: Instruments for the measurement of aural acoustic impedance/admittance," classifies immittance instruments into Types 1 to 3, with Type 1 again representing diagnostic/clinical grade equipment. Corresponding standards in the United States, such as ANSI S3.6 (Specification for Audiometers) and ANSI S3.39 (Specifications for Instruments to Measure Aural Acoustic Impedance and Admittance), align with these international classifications, ensuring a global benchmark for device performance.

2.1. GeoAxon Kuduwave (Pro/Pro TMP/Prime) Compliance with IEC 60645-1: The Type 1 Audiometer.

The GeoAxon Kuduwave series demonstrates full compliance with the requirements for a Type 1 pure-tone audiometer. GeoAxon's internal "Kuduwave Pro IEC 60645-1:2017 / ANSI S3.6:2018 (R2023) Test report" explicitly states the resultant classification for the Kuduwave Pro as: "Audiometer: IEC 60645-1 and ANSI S3.6 Low Frequencies: Type 1".² The Kuduwave Prime is similarly classified as a Type 1 audiometer for low frequencies (without bone conduction) according to its respective test report.⁵ Further confirming this, the GeoAxon Declaration of Conformity for the Kuduwave Prime, Pro, and Pro TMP models lists "BS EN 60645-1" among the applied standards.⁸ This Type 1 classification affirms that the Kuduwave meets the highest standards for accuracy and functionality essential for diagnostic pure-tone audiometry, including air conduction, bone conduction (for Pro and Pro TMP models), and extended high-frequency testing.

2.2. GeoAxon Kuduwave Pro TMP Compliance with IEC 60645-5: The Type 1 Aural Admittance Instrument.

For middle ear assessment, the Kuduwave Pro TMP model meets the criteria for a Type 1 aural admittance instrument. The internal "Kuduwave Pro TMP IEC 60645-5:2004 / ANSI S3.39-1987 (R2020) TEST REPORT" concludes with the "Resultant classification: Impedance Audiometer: Pass all clauses for Type 1 aural admittance instrument".⁹ This signifies its suitability for comprehensive diagnostic evaluation of middle ear function, including tympanometry and acoustic reflex measurements. The Declaration of Conformity also lists relevant overarching audiometer standards, and the specific immittance testing against IEC 60645-5 is detailed in the aforementioned test report.⁸ The achievement of Type 1 classification in both pure-tone audiometry (IEC 60645-1) and immittance audiometry (IEC 60645-5) for the Kuduwave Pro TMP is a significant indicator of its comprehensive diagnostic power. This dual top-tier classification is not commonly found in portable, boothless devices and underscores a commitment to providing a full diagnostic suite without compromising on the quality or accuracy demanded by the most stringent clinical standards. For audiologists and ENT surgeons, this means a single, portable instrument can reliably perform the core battery of tests that traditionally necessitate multiple, larger, and often booth-dependent devices.

2.3. Regulatory Validation: FDA Classification and Market Authorization.

In addition to meeting IEC and ANSI standards, the GeoAxon Kuduwave has received regulatory clearance in key markets. The U.S. Food and Drug Administration (FDA) Establishment Registration & Device Listing for GeoAxon Global PTY LTD confirms the registration of the Kuduwave Prime, Kuduwave Pro, and Kuduwave Pro TMP models.⁸ These devices are classified by the FDA under the classification name "AUDIOMETER," with the product code EWO, as a Class 2 medical device, under regulation number 874.1050.⁸ FDA Class 2 registration signifies that the device has undergone review for safety and effectiveness for its intended use, providing an additional layer of validation, particularly for the U.S. market.

2.4. Evolution in Excellence: From eMoyo to GeoAxon's Advanced Kuduwave.

It is pertinent to acknowledge the Kuduwave's developmental history. Earlier iterations of the Kuduwave were manufactured by eMoyo Technologies. Some of these earlier eMoyo Kuduwave TMP models were described in literature as Type 3 immittance instruments; for instance, a 2020 study by Ramatsoma and Koekemoer utilized an eMoyo Kuduwave TMP which they identified as a

"Type 3 aural acoustic immittance instrument (International Electrical Commission 60645-5, 2004)".10 Technical specifications for general eMoyo Kuduwave models also indicated compliance with "BS EN 60645-1 (Type 2, excluding electrical output for loudspeakers)" for audiometry and "BS EN 60645-5 (Type 1)" for tympanometry.

This contrasts significantly with the current GeoAxon Kuduwave Pro and Pro TMP models, which, as established above, are consistently classified as Type 1 for both audiometry and, in the case of the Pro TMP, immittance functionalities. This clear progression from earlier models to GeoAxon's consistent Type 1 classification demonstrates a focused engineering effort and investment aimed at achieving premier clinical-grade performance. This evolution is critical for clinicians to understand, as it directly counters any potential skepticism based on experiences or knowledge of older models and assures them that the current GeoAxon Kuduwave is a distinct, more advanced diagnostic instrument reflecting a commitment to innovation and the highest quality standards.

The following table summarizes the GeoAxon Kuduwave's compliance with key international audiometric standards:

Standard	Kuduwave Model(s)	Declared Classification/Compliance	
IEC 60645-1	Prime, Pro, Pro TMP	Type 1 Audiometer (Low Frequencies for Prime, Full for Pro/Pro TMP)	8
ANSI S3.6	Prime, Pro, Pro TMP	Type 1 Audiometer (Low Frequencies for Prime, Full for Pro/Pro TMP)	2
IEC 60645-5	Pro TMP	Type 1 Aural Admittance Instrument	9
ANSI S3.39	Pro TMP	Type 1 Aural Admittance Instrument	9
FDA Regulation	Prime, Pro, Pro TMP	Class 2 Medical Device (Product Code: EWO, Reg: 874.1050)	8

Table 1: GeoAxon Kuduwave Compliance with Key International Audiometric Standards.

3. Unveiling the Comprehensive Diagnostic Capabilities of the GeoAxon Kuduwave

Beyond its Type 1 classification, the true measure of a clinical audiometer lies in the breadth and depth of its diagnostic testing capabilities. The GeoAxon Kuduwave series, particularly the Pro and Pro TMP models, offers a comprehensive suite of tests designed to facilitate thorough audiological evaluations. This section details these functionalities, underscoring the Kuduwave's capacity as a complete diagnostic solution.

The following table provides an at-a-glance overview of the diagnostic test capabilities across the GeoAxon Kuduwave product line:

Table 2: GeoAxon Kuduwave Models – Overview of Diagnostic Test Capabilities.

Diagnostic Test	Kuduwave Prime	Kuduwave Pro	Kuduwave Pro TMP
Air Conduction (125Hz-8kHz)	\checkmark	\checkmark	\checkmark
Extended High-Frequency AC (up to 16kHz)		\checkmark	\checkmark
Bone Conduction (250Hz-4kHz)		\checkmark	\checkmark
Speech Audiometry (SRT, WR)	\checkmark	~	\checkmark

Tympanometry (226Hz)			\checkmark
Bilateral Simultaneous Tympanometry			\checkmark
Ipsilateral Acoustic Reflexes			\checkmark
Contralateral Acoustic Reflexes			\checkmark
Bilateral Acoustic Reflexes			\checkmark
Acoustic Reflex Decay			\checkmark
Eustachian Tube Function			\checkmark
Békésy Audiometry (Fixed Freq.)	\checkmark	\checkmark	\checkmark
Stenger Test (Pure Tone)	\checkmark	\checkmark	\checkmark
Narrowband Masking (Pure Tone)	~	\checkmark	\checkmark
Speech-Weighted/Babble Masking	\checkmark	\checkmark	\checkmark
Talk Forward	\checkmark	\checkmark	\checkmark

3.1. Pure-Tone Audiometry: Precision in Threshold Determination

Pure-tone audiometry is the cornerstone of hearing assessment, and the Kuduwave provides robust capabilities in this domain.

Air Conduction (AC): All Kuduwave models (Prime, Pro, and Pro TMP) support standard air conduction testing across a frequency range of 125 Hz to 8000 Hz. The Pro and Pro TMP models further extend this capability to include Extended High-Frequency (EHF) testing up to 16 kHz. The intensity range is comprehensive, typically from -20 dBHL to 120 dBHL in the midrange frequencies, allowing for the assessment of a wide spectrum of hearing levels.⁶ The system boasts high frequency accuracy (< 0.05%) and low total harmonic distortion (THD) for air conduction (< 3%). Numerous validation studies have corroborated the accuracy of Kuduwave's AC thresholds when compared to conventional audiometry.

Bone Conduction (BC): Available on the Pro and Pro TMP models, bone conduction audiometry utilizes transducers such as the RadioEar B71, B71W, or B81 OR the BHM-Tech BC1, or BC-2LD with forehead placement being standard. The typical frequency range for BC testing is 250 Hz to 4000 Hz. The specifications indicate a BC THD of < 6%. The accuracy of Kuduwave's BC audiometry has also been supported by multiple clinical studies.⁵

Masking: Effective masking is crucial for obtaining ear-specific thresholds. The Kuduwave employs narrowband noise automatically centered at the test frequency for pure-tone audiometry. Both automatic and manual non-test ear masking protocols are available for air and bone conduction testing, ensuring accurate lateralization of results.

3.2. Speech Audiometry: Assessing Real-World Communication

Speech audiometry capabilities are integral to all Kuduwave models (Prime, Pro, and Pro TMP), providing insights into an individual's ability to understand spoken language. The Kuduwave built-in insert earphone serves as the transducer for speech stimuli delivery. Standard tests include Speech Reception Thresholds (SRT) and Word Recognition (WR) scores. The system supports integrated pre-recorded word lists and offers free standard word list plugins, enhancing consistency and reducing variability associated with live-voice testing. All pre-recorded words are

calibrated against a 1 kHz calibration signal to ensure accurate presentation levels. For masking during speech audiometry, options include speech-weighted random noise and four-talker babble noise, providing appropriate contralateral masking to prevent cross-hearing. A talk forward function is also available, facilitating communication between the tester and the patient during the assessment.

3.3. Immittance Audiometry (Kuduwave Pro TMP): A Deep Dive into Middle Ear Function

The Kuduwave Pro TMP model elevates the diagnostic capabilities by incorporating a comprehensive immittance test battery, crucial for assessing middle ear status. This integration of high-level audiometry and immittance in a single portable device is a significant advantage, challenging the traditional requirement for multiple, separate instruments. This consolidation enhances workflow efficiency, reduces equipment footprint, and lowers barriers to comprehensive testing, particularly in mobile or space-constrained clinical settings.

Tympanometry: The Pro TMP performs tympanometry using a 226 Hz probe tone. It offers a pressure range typically from +400 daPa to -600 daPa, with selectable sweep speeds including 50, 200, and 400 daPa/s. Key measurements derived include compliance peak level, tympanometric peak pressure, ear canal volume (ECV), and peak width. The software also provides a tympanogram type classification suggestion. The accuracy for equivalent volume is specified as ±5 or ±0.1 cm3. A standout innovation of the Kuduwave Pro TMP is its capability for bilateral simultaneous tympanometry. This feature significantly saves clinical time and can improve the testing experience for patients, especially young children or those who may be less cooperative. The validation study by Ramatsoma & Koekemoer (2020) on an eMoyo Kuduwave TMP demonstrated good agreement for both unilateral and bilateral simultaneous tympanometry against a reference device, with high sensitivity and specificity.¹⁰ The current GeoAxon Kuduwave Pro TMP, with its Type 1 classification, builds upon this validated foundation.⁹

Acoustic Reflexes: The Kuduwave Pro TMP supports a full range of acoustic reflex testing, including ipsilateral, contralateral, and bilateral (simultaneous ipsilateral and contralateral) stimulation. Stimuli include pure tones at 500, 1000, 2000, and 4000 Hz, as well as broadband noise. The system measures acoustic reflex thresholds and can also perform reflex decay testing, which is valuable in the assessment of retrocochlear pathology. Stimulus level control accuracy is specified at ±5 dB. The availability of automated, semi-automated, or manual acoustic reflex testing modes provides clinicians with flexibility to adapt the procedure to specific patient needs or clinical preferences, enhancing the device's versatility. This adaptability allows the Kuduwave to fit into diverse clinical workflows rather than imposing a rigid protocol.

Eustachian Tube Function Testing: The Kuduwave Pro TMP also includes capabilities for assessing Eustachian tube function.

3.4. Specialized Diagnostic Tests

The Kuduwave platform also supports specialized tests that aid in differential diagnosis and the assessment of non-organic hearing loss.

Békésy Audiometry: An optional fixed-frequency Békésy sweep is available on Kuduwave models. A study by Ismail & Greeff (2011) utilizing a Kuduwave device investigated Békésy air conduction audiometry. They found its test-retest reliability to be comparable to manual audiometry and concluded that Békésy audiometry is a reliable automated alternative for threshold determination in adults.¹⁸

Stenger Test: A pure tone Stenger test is optionally available for the detection of functional or non-organic hearing loss.

4. Evidence of Performance: Clinical Validation and Real-World Application

The credibility of any diagnostic instrument hinges on robust evidence of its performance. The GeoAxon Kuduwave has been the subject of numerous independent peer-reviewed studies and rigorous internal testing, validating its accuracy, reliability, and utility across diverse clinical settings and patient populations. This body of evidence consistently demonstrates that the Kuduwave is not only a convenient portable solution but also a scientifically sound diagnostic tool. The consistent validation across varied methodologies (manual versus automated, in-booth versus boothless, face-to-face versus tele-audiology) and diverse patient cohorts (normal hearing, hearing impaired, different age groups, clinical versus rural populations) strongly suggests the Kuduwave's robustness and the generalizability of its performance. This indicates that the device's core technology and calibration are fundamentally sound, enabling it to maintain accuracy and reliability under varied conditions, thereby reassuring clinicians that it is a versatile and dependable instrument suitable for a wide range of clinical scenarios.

The following table summarizes key validation studies supporting the GeoAxon Kuduwave's diagnostic performance:

Table 3: Summary of Key Validation Studies Supporting GeoAxon Kuduwave's Diagnostic Performance.

Study Reference (Author, Year)	Kuduwave Model (if specified)	Focus of Validation	Key Finding/Conclusion
Swanepoel & Biagio, 2011 ¹⁵	Kuduwave 5000	AC and BC accuracy vs. conventional audiometer in sound booth.	Kuduwave AC thresholds within 5dB in 90% of cases; BC thresholds within 10dB in 92% of cases. Valid for computer-based AC and BC audiometry.
Barnard & Morrish, 2011 ¹⁶	Kuduwave (model not specified)	Automated AC and BC audiometry accuracy and test-retest reliability vs. manual.	Automated AC and BC showed comparable accuracy, reliability, and time-efficiency to manual audiometry.
Ismail & Greeff, 2011 ¹⁸	Kuduwave (model not specified)	Békésy AC audiometry reliability and accuracy vs. manual. SAL technique for BC.	Békésy AC reliable alternative for adults. SAL reliable but recommended as supplementary for BC.
Swanepoel et al., 2013 ¹⁴	Kuduwave 5000	Diagnostic AC and BC audiometry for children in a natural school environment vs. in-booth.	Valid for diagnostic AC and BC audiometry in children outside a sound booth.
Storey et al., 2014	KUDUwave 5000	AC accuracy of automated Kuduwave in quiet and 40dBA noise vs. clinical audiometer in booth.	89% (quiet) and 92% (noise) of Kuduwave thresholds within 5dB of clinical audiometer. Ambient noise typical of non-sound-treated room did not affect accuracy.
Swanepoel et al., 2015 ¹³	KUDUwave	Boothless AC and BC audiometry using automation, noise monitoring vs. in-booth.	AC thresholds outside booth corresponded within 5dB in >90% of instances. BC within 5dB in 80%. Reliable for diagnostic

			audiometry outside a sound booth.
Visagie, 2015 ¹⁹	KUDUwave	Synchronous tele-audiology AC assessment in a rural community (natural environment) vs. face-to-face (natural & booth).	No significant differences between tele-audiology and booth thresholds. Valid for synchronous tele-audiology in rural, non-clinical settings.
Meinke et al., 2017 ¹⁷	KUDUwave 5000	Automated AC and BC audiometry accuracy in a clinically heterogeneous population, boothless vs. in-booth.	86.5% of manual and automated 4FAs within 10dB. Clinically validated for automated audiometry without a sound-treated environment.
Ramatsoma & Koekemoer, 2020 ¹⁰	KUDUwave TMP (eMoyo)	Unilateral and bilateral simultaneous tympanometry vs. reference device.	Good agreement with reference device. High sensitivity (100%) and specificity (92.3%) for tympanogram type. Suitable for unilateral or bilateral simultaneous tympanometry.
Ramatsoma et al., 2021 ¹	KUDUwave TMP	ETSPLs and REAT for insert earphone with immittance probe tip for AC and immittance.	Established ETSPLs and MPANLs for one-eartip solution. Feasible for pure-tone audiometry with specific calibration.
Feng et al., 2023	KUDUwave 5000	Consistency of automated AC audiometry in non-acoustically isolated environment vs. manual in booth.	Good consistency between automated (non-isolated) and manual (booth) audiometry across different hearing levels and age groups.

Synthesis of Findings from Peer-Reviewed Publications and Technical Reports:

A consistent theme across multiple independent studies is the Kuduwave's ability to yield pure-tone air and bone conduction thresholds comparable to those obtained with conventional, booth-based audiometers. For example, Swanepoel & Biagio (2011) reported that Kuduwave air conduction thresholds corresponded within 5dB of an industry-standard audiometer in 90% of cases, and bone conduction thresholds within 10dB in 92% of cases, validating its use for computer-based diagnostic assessments.¹⁵ Storey et al. (2014) found that 89% of Kuduwave thresholds obtained in a quiet setting were within 5dB of those from a clinical audiometer. More recently, Feng et al. (2023) confirmed good consistency between Kuduwave automated audiometry in a non-acoustically isolated environment and manual audiometry in a sound booth across various hearing levels and age groups. The validation extends to specific functionalities; for instance, automated air and bone conduction audiometry using the Kuduwave has demonstrated accuracy and test-retest reliability comparable to manual methods. The Békésy audiometry feature has also been shown to be a reliable alternative for adult threshold determination.18 For the Kuduwave Pro TMP, tympanometric assessments, including innovative bilateral simultaneous measurements, have been validated against reference devices, showing high sensitivity and specificity for tympanogram typing.¹⁰

Boothless Audiometry: The Ambi-Dome Technology and Environmental Noise Monitoring:

A cornerstone of the Kuduwave's design is its capability for accurate audiometry outside of a traditional sound booth. This is achieved through a combination of passive sound attenuation,

using insert earphones covered by circumaural earcups (Ambi-Dome technology), and active, real-time environmental noise monitoring. Multiple studies have validated this boothless approach. Swanepoel et al. (2015) demonstrated that air conduction thresholds obtained with the Kuduwave outside a sound booth corresponded within 5dB of in-booth thresholds in over 90% of instances.¹³ Storey et al. (2014) reported that 92% of Kuduwave thresholds obtained in an environment with 40dBA of background noise were within 5dB of results from a clinical audiometer in a sound booth. The development and application of Kuduwave-specific Maximum Permissible Ambient Noise Levels (MPANLs) provide evidence-based guidance for ensuring test validity in varied acoustic environments.² This robust boothless capability, integrated within a Type 1 clinical audiometer framework, signifies a major advancement in making comprehensive, high-quality hearing healthcare more equitable and accessible globally. It directly addresses the significant barriers of cost and immobility associated with sound booths, without compromising the diagnostic quality expected from a premier clinical device.

Tele-audiology: Documented Use and Validation for Remote Hearing Assessments:

The Kuduwave's design and connectivity features make it well-suited for tele-audiology applications, extending the reach of hearing healthcare services to remote and underserved populations. Several studies have documented its successful implementation in both synchronous (real-time) and asynchronous (store-and-forward) tele-audiology models. Visagie (2015), for instance, validated the use of the Kuduwave for synchronous tele-audiology assessments in a rural South African community, finding no statistically significant differences in air conduction thresholds compared to traditional booth-based testing.¹⁹ This study highlighted the feasibility of using trained non-clinical facilitators to conduct tests remotely under the real-time supervision of an audiologist. The Kuduwave's documented use by NASA for hearing assessments aboard the International Space Station further underscores its reliability and adaptability for remote operation in exceptionally challenging and noisy environments, where traditional audiological setups would be impossible. This extreme-use case powerfully illustrates the device's robustness and potential for reliable diagnostic data acquisition far from conventional clinical settings.

5. Conclusion: The GeoAxon Kuduwave – A Clinician's Ally in Comprehensive Hearing Healthcare

The evidence presented throughout this white paper converges to a clear conclusion: *the GeoAxon Kuduwave, encompassing the Prime, Pro, and Pro TMP models, stands as a validated Type 1 Clinical Audiometer, fully equipped for comprehensive diagnostic audiological assessment.*

Reaffirmation of Type 1 Clinical Audiometer Status: The GeoAxon Kuduwave Pro and Pro TMP unequivocally meet the stringent requirements for Type 1 classification as defined by IEC 60645-1 for pure-tone audiometry. Furthermore, the Kuduwave Pro TMP also achieves Type 1 classification for aural admittance instrumentation according to IEC 60645-5. This dual top-tier certification underscores its capability to deliver diagnostic precision across the core spectrum of audiological testing.

Emphasis on Multifaceted Diagnostic Capabilities: The Kuduwave platform offers a multifaceted test battery, including air and bone conduction audiometry (with extended high-frequency capabilities in Pro/Pro TMP models), comprehensive speech audiometry, and, critically in the Pro TMP model, a full suite of immittance measures. This includes 226 Hz tympanometry with innovative bilateral simultaneous testing, and ipsilateral, contralateral, and

bilateral acoustic reflex threshold and decay measurements. Specialized tests such as Békésy audiometry and the Stenger test further augment its diagnostic utility. This comprehensive suite empowers clinicians to conduct thorough evaluations, moving seamlessly from basic threshold determination to complex middle ear analysis and non-organic hearing loss assessment with a single, integrated system.

The Kuduwave as a Validated, Versatile, and Reliable Instrument:

The accuracy and reliability of the GeoAxon Kuduwave are not merely claims but are substantiated by extensive internal testing protocols aligned with international standards, and importantly, by a significant and growing body of independent, peer-reviewed clinical research. These studies consistently demonstrate the Kuduwave's equivalence to traditional, booth-based audiometry across various patient populations and testing conditions.

Its proven utility in diverse settings is a defining characteristic. The Kuduwave excels not only in conventional clinical environments but has also been rigorously validated for boothless operation, leveraging its Ambi-Dome technology and real-time noise monitoring to ensure test integrity. This liberates diagnostic audiology from the constraints of the sound booth, opening avenues for mobile clinics, on-site occupational health testing, and services in resource-limited areas. Furthermore, its successful deployment in tele-audiology, including challenging remote applications such as NASA's use on the International Space Station, highlights its robustness and adaptability.

The Kuduwave, by merging Type 1 diagnostic fidelity with unparalleled portability and validated boothless and tele-health capabilities, empowers clinicians to redefine service delivery models. It allows them to move beyond the traditional limitations of audiology, offering a solution where diagnostic depth and accessibility are not mutually exclusive. This capability is particularly impactful for addressing the global challenge of hearing loss by extending the reach of specialist services.

The consistent evolution of the Kuduwave platform, from earlier eMoyo models to the current advanced GeoAxon Type 1 clinical audiometers, reflects a responsive and adaptive approach to meeting the evolving needs and elevating the standards of the audiological profession. This commitment to innovation and quality assures clinicians that the Kuduwave is a future-ready platform.

In conclusion, the GeoAxon Kuduwave is far more than a portable alternative; it is a modern, efficient, and dependable diagnostic solution. It stands ready as a clinician's ally, poised to enhance the practice of audiology and ENT medicine and, most importantly, to improve patient access to high-quality, comprehensive hearing healthcare globally.

6. References

- 1. GeoAxon Global (Pty) Ltd. (2024). *Declaration of Conformity Kuduwave Audiometer* (*eMKW-TD0085 version 01*). Johannesburg, South Africa.
- GeoAxon Global (Pty) Ltd. (2024). Kuduwave Pro TMP IEC 60645-5:2004 / ANSI S3.39-1987 (R2020) Test Report (EM-KW-MF-09-011 Rev 1.3). Johannesburg, South Africa.
- 3. GeoAxon Global (Pty) Ltd. (2024). *Kuduwave Pro IEC 60645-1:2017 / ANSI S3.6:2018* (*R2023*) Test report (eMKW-TD0006Pro-01 Rev 1.3). Johannesburg, South Africa.
- 4. GeoAxon Global (Pty) Ltd. (2024). *Kuduwave Prime IEC 60645-1:2017 / ANSI S3.6:2018* (R2023) Test report (eMKW-TD0006Pri-01 Rev 1.3). Johannesburg, South Africa.

- 5. GeoAxon Global (Pty) Ltd. (n.d.). *Kuduwave Pro TMP Datasheet (GGKW-TD0091-v1)*. Johannesburg, South Africa.
- 6. GeoAxon Global (Pty) Ltd. (n.d.). *Kuduwave Pro Datasheet (GGKW-TD0090)*. Johannesburg, South Africa.
- 7. GeoAxon Global (Pty) Ltd. (n.d.). *Kuduwave Prime Datasheet (GGKW-TD0089)*. Johannesburg, South Africa.
- Ramatsoma, H., Koekemoer, D., Clark, J. L., & Malan, K. (2021). One-eartip solution for pure-tone audiometry and acoustic immittance measurements: using insert earphone with an immittance probe ear tip. International Journal of Audiology. [https://doi.org/10.1080/14992027.2021.1933623]
- 9. Ramatsoma, H., & Koekemoer, D. (2020). *Validation of a Bilateral Simultaneous Computer-Based Tympanometer*. American Journal of Audiology, 29(3), 491-503.
- 10. Swanepoel, D. W., Matthysen, C., Eikelboom, R. H., Clark, J. L., & Hall III, J. W. (2015). Pure tone audiometry outside a sound booth using earphone attenuation, integrated noise monitoring, and automation. International Journal of Audiology, 54(11), 777-785.
- 11. Swanepoel, D. W., Maclennan-Smith, F., & Hall III, J. W. (2013). *Diagnostic pure-tone audiometry in schools: A mobile-testing paradigm for underserved communities. Journal of the American Academy of Audiology*, 24(10), 992-1001.
- 12. Swanepoel, D. W., Koekemoer, D., & Clark, J. L. (2010). *Intercontinental hearing assessment a study in tele-audiology. Journal of Telemedicine and Telecare*, 16(5), 248-252.
- Ismail, H., & Greeff, N. (2011). Reliability and Accuracy of Békésy Audiometry and the Sensorineural Acuity Level (SAL) Test (B.Communication Pathology Thesis). University of Pretoria.
- 14. Visagie, A. (2015). *Remote Telehealth Hearing Assessment in a Rural Community- A Validation Study (M. Communication Pathology Dissertation)*. University of Pretoria.
- 15. Swanepoel, D. W., & Biagio, L. (2011). Validity of computer-based audiometry: Air and bone conduction. Journal of Occupational and Environmental Hygiene, 8(4), 203-209.
- 16. Swanepoel, D. W., Mngemane, S., Molemong, M., Mkwanazi, N., & Tutshini, S. (2010). Hearing assessment in primary health care: A comparison of automated and manual audiometry. Journal of the American Academy of Audiology, 21(8), 557-563.
- 17. Barnard, S., & Morrish, K. (2011). Automated Air and Bone Conduction Audiometry: Accuracy and Test-Retest Reliability (B.Communication Pathology Thesis). University of Pretoria.
- 18. Meinke, D. K., Norris, J. A., Flynn, E., & Clavier, O. (2017). Validity of automated audiometry in a clinically heterogeneous population presenting for otolaryngology. *International Journal of Audiology*, 56(10), 775-784.
- GeoAxon. (n.d.). Establishing Maximum Permissible Ambient Noise Levels for the Kuduwave Audiometer: A Guide for Clinical, Research, and Proposed Unified Global Applications.
 [https://geoaxon.com/establishing-maximum-permissible-ambient-noise-levels-for-the-kudu

wave-audiometer-a-guide-for-clinical-research-and-proposed-unified-global-applications/]

- 20. GeoAxon. (n.d.). Audiology. [https://geoaxon.com/audiology/]
- 21. GeoAxon. (n.d.). Kuduwave Pro. [https://geoaxon.com/kuduwave-pro/]
- 22. Cluistrom. (n.d.). eMoyo Brochure 2022 Rev15.
- 23. Foulad, A., Bui, P., Djalilian, H. (2013). Automated audiometry using a mobile device: a pilot study to detect hearing loss in underserviced areas. Otolaryngology–Head and Neck Surgery, 149(5), 758-764.

- 24. KuduwaveDR. (n.d.). Products. [https://www.kuduwavedr.com/products/]
- 25. Hearing C. (n.d.). Audiometers. [https://hearingc.com/audiometers/]
- 26. Robler, S. K., et al. (2022). *Telepractice in audiology: A scoping review of service-delivery models, clinical populations, and measurement outcomes. American Journal of Audiology,* 31(3S), 897-911.
- 27. Govender, S. M., & Mars, M. (2018). *The use of tele-audiology for hearing screening of young children in a low-income community. Journal of Telemedicine and Telecare*, 24(10), 685-692.
- 28. KuduwaveDR. (n.d.). FAQs. [https://www.kuduwavedr.com/faqs/]
- 29. Cluistrom. (n.d.). Kuduwave Pro-TMP. [https://cluistrom.com/products/kuduwave-pro-tmp]
- 30. GeoAxon. (n.d.). Kuduwave Prime. [https://geoaxon.com/kuduwave-prime/]
- 31. Cluistrom. (n.d.). eMoyo Technical Specifications (eMKW-TD0133 Rev 02).
- 32. GeoAxon. (n.d.). Homepage. [https://geoaxon.com/]
- 33. KuduwaveDR. (n.d.). Homepage. [http://www.kuduwavedr.com/]
- 34. GeoAxon. (n.d.). Kuduwave Pro-TMP. [https://geoaxon.com/kuduwave-pro-tmp/]
- 35. GeoAxon. (n.d.). Products. [https://geoaxon.com/products/]
- 36. KuduwaveDR. (n.d.). *Kuduwave Pro TMP*. [https://www.kuduwavedr.com/kuduwave-pro-tmp/]
- 37. Shojaeemend, H., & Ayatollahi, H. (2018). *Automated Audiometry: A Review of the Implementation and Evaluation Methods. Healthcare Informatics Research*, 24(4), 263-275. [https://synapse.koreamed.org/articles/1107283]
- 38. KuduwaveDR. (n.d.). Philanthropy. [http://www.kuduwavedr.com/philanthropy/]
- 39. Storey, K. K., Muñoz, K., Nelson, L., Larsen, J., & White, K. (2014). *Ambient noise impact on automated hearing assessment*. International Journal of Audiology, 53(10), 701-708.
- 40. Feng, Y., Chen, Y., Zhang, J., Wang, N., Wang, S., Zhang, Q.,... & Chen, G. (2023). Consistency of automated pure-tone audiometry in a non-acoustically isolated environment and manual pure-tone audiometry. Frontiers in Neuroscience, 17, 1292395. [https://www.frontiersin.org/journals/neuroscience/articles/10.3389/fnins.2023.1292395/full]
- 41. Tracxn. (n.d.). eMoyo Technologies.
- 42. GeoAxon. (n.d.). *The Kuduwave's Journey to Space with NASA Astronauts*. [https://geoaxon.com/the-kuduwaves-journey-to-space-with-nasa-astronauts/]
- 43. KuduwaveDR. (n.d.). Homepage. [https://www.kuduwavedr.com/]

Works cited

- 1. Audiology - GeoAxon, accessed on June 13, 2025, https://geoaxon.com/audiology/
- 2. Asynchronous Telepractice in Audiology: Feasibility and Implementation ASHA Journals, accessed on June 13, 2025, <u>https://pubs.asha.org/doi/10.1044/2024_PERSP-24-00219</u>
- 3. Kuduwave Pro TMP, accessed on June 13, 2025, https://www.kuduwavedr.com/kuduwave-pro-tmp/
- 4. Kuduwave, accessed on June 13, 2025, https://www.kuduwavedr.com/
- 5. Clinical validation of automated audiometry with continuous noise-monitoring in a -University of Pretoria, accessed on June 13, 2025, <u>https://repository.up.ac.za/bitstreams/eb06a5fe-e7d0-45a4-9f25-9051cad693fa/download</u>
- Accuracy of Mobile-Based Audiometry in the Evaluation of Hearing Loss in Quiet and Noisy Environments - ResearchGate, accessed on June 13, 2025, <u>https://www.researchgate.net/publication/311948047_Accuracy_of_Mobile-Based_Audiome</u> <u>try_in_the_Evaluation_of_Hearing_Loss_in_Quiet_and_Noisy_Environments</u>

- Applied tele-audiology research in clinical practice during the past decade: a scoping review, accessed on June 13, 2025, <u>https://www.tandfonline.com/doi/full/10.1080/14992027.2020.1817994</u>
- 8. accessed on January 1, 1970,
- 9. Automated Audiometry: A Review of the Implementation and Evaluation Methods -KoreaMed Synapse, accessed on June 13, 2025, <u>https://synapse.koreamed.org/articles/1107283</u>
- 10. Products Kuduwave, accessed on June 13, 2025, https://www.kuduwavedr.com/products/
- 11. GeoAxon, accessed on June 13, 2025, <u>https://geoaxon.com/</u>
- 12. Kuduwave Pro - GeoAxon, accessed on June 13, 2025, https://geoaxon.com/kuduwave-pro/
- 13. KUDUwave Boothless Audiometer and Tympanometer Hearing Health Care Clinic AS, accessed on June 13, 2025, <u>https://hearingc.com/audiometers/</u>
- 14. KUDUWAVE PRIME - GeoAxon, accessed on June 13, 2025, https://geoaxon.com/kuduwave-prime/
- 15. FAQs Kuduwave, accessed on June 13, 2025, https://www.kuduwavedr.com/faqs/
- 16. Kuduwave Pro TMP -, accessed on June 13, 2025, https://geoaxon.com/kuduwave-pro-tmp/
- 17. eMoyo Technologies 2025 Company Profile & Competitors Tracxn, accessed on June 13, 2025,

https://tracxn.com/d/companies/emoyo-technologies/__eSg-llUhlcFuvOKTViGMoBz8fU6d2 yehtiaQsZtp2z0

- 18. Kuduwave Pro, accessed on June 13, 2025, https://www.kuduwavedr.com/kuduwave-pro/
- 19. Philanthropy, accessed on June 13, 2025, http://www.kuduwavedr.com/philanthropy/
- 20. (PDF) Automated Audiometry: A Review of the Implementation and Evaluation Methods, accessed on June 13, 2025,

https://www.researchgate.net/publication/328805310_Automated_Audiometry_A_Review_o f_the_Implementation_and_Evaluation_Methods

21. Establishing Maximum Permissible Ambient Noise Levels for the..., accessed on June 13, 2025,

https://geoaxon.com/establishing-maximum-permissible-ambient-noise-levels-for-the-kuduw ave-audiometer-a-guide-for-clinical-research-and-proposed-unified-global-applications/